

## Product Datasheet

### BI-1347 (orb573443)

<b>Description</b>	BI-1347 is a potent, selective inhibitor of CDK8/cyclinC with IC50 of 1 nM.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{22}H_{20}N_4O$
<b>MW</b>	356.429
<b>CAS Number</b>	[2163056-91-3]
<b>SMILES</b>	<chem>CN(C)C(=O)CN1C=C(C=N1)C2=CC=C(C=C2)C3=CN=CC4=CC=CC=C43</chem>
<b>Expiration Date</b>	12 months from date of receipt.