



Product Datasheet

HSP90AA1 antibody (orb573596)


Description	Rabbit polyclonal antibody to HSP90AA1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human HSP90AA1
Target	HSP90AA1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: HLYKDLQPFILLRLLMPEETQTQDQPMEEVEVETFAFQAEIAQLMSLIIN
Clonality	Polyclonal
MW	98kDa
Uniprot ID	P07900
NCBI	NP_001017963
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
HSP90AA1,
Sam...



Anti-HSP90AA1
/ Hsp90
antibody IHC
stain...



Anti-HSP90AA1
/ Hsp90
antibody IHC
stain...