



Product Datasheet

Nkx2-5 antibody (orb573852)


Description	Rabbit polyclonal antibody to Nkx2-5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the n terminal region of human Nkx2-5
Target	Nkx2-5
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ELGRAPSPAKCASAFPAAPAFYPRAYSDDPAKDPRAEKKELCALQKAVE
Clonality	Polyclonal
MW	35kDa
Uniprot ID	P52952
NCBI	NP_004378
Expiration Date	12 months from date of receipt.



Host: Mouse,
Target Name:
NKX2-5,
Sample...



Rabbit Anti-
NKX2-5
Antibody,
Catalog Num...



WB Suggested
Anti-Nkx2-5
Antibody
Titrat...