



Product Datasheet

CSEN antibody (orb573970)

Description	Rabbit polyclonal antibody to CSEN
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human CSEN
Target	KCNIP3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: STVRHQPEGLDQLQAQTKFTKKELQSLYRGFKNECPTGLVDEDTFKLIYA
Clonality	Polyclonal
MW	29kDa
Uniprot ID	Q9Y2W7
NCBI	NP_038462
Expiration Date	12 months from date of receipt.



Host: Mouse,
Target Name:
KCNIP3,
Sample...



Human kidney



WB Suggested
Anti-CSEN
Antibody
Titratio...