



Product Datasheet

SIRT5 antibody (orb574179)


Description	Rabbit polyclonal antibody to SIRT5
Species/Host	Rabbit
Reactivity	Bovine, Goat, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human SIRT5
Target	SIRT5
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LSGKGAPEPGTQDASIPVEKLPRCEEAGCGLLRPHVVWFGENLDPAIL
Clonality	Polyclonal
MW	33kDa
Uniprot ID	Q9NXA8
NCBI	NP_112534
Expiration Date	12 months from date of receipt.



Human kidney



Human Lung



Rabbit Anti-SIRT5 Antibody, Paraffin Emb...