



Product Datasheet

NCOA3 antibody (orb574233)

Description	Rabbit polyclonal antibody to NCOA3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human NCOA3
Target	NCOA3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DQNPVESSMCQSNRDLSDKESKESSEVEGAENQRGPLESKGGHKKLLQLL
Clonality	Polyclonal
MW	155kDa
Uniprot ID	Q0VF45
NCBI	NP_006525
Expiration Date	12 months from date of receipt.



Host: Rabbit, Target Name: NCOA3, Sample...



Immunohistochemistry with Human Tonsil I...



WB Suggested Anti-NCOA3 Antibody Titrati...