



Product Datasheet

PCYOX1 antibody (orb574476)

Description	Rabbit polyclonal antibody to PCYOX1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human PCYOX1
Target	PCYOX1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: IFSQETLTKAQILKFLSYDYAVKKPWLAYPHYKPPKPCPSIILHDLRY
Clonality	Polyclonal
MW	57kDa
Uniprot ID	Q9UHG3
NCBI	NP_057381
Expiration Date	12 months from date of receipt.



Human Liver



Rabbit Anti-PCYOX1 Antibody, Paraffin Em...



WB Suggested Anti-PCYOX1 Antibody Titrat...