



Product Datasheet

PLD2 antibody (orb574485)


Description	Rabbit polyclonal antibody to PLD2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human PLD2
Target	PLD2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DLHYRLTDLGDSSESAASQPPTPRPDPSPATPDLSHNQFFWLKGDYSNLIT
Clonality	Polyclonal
MW	106 kDa
Uniprot ID	O14939
NCBI	NP_002654
Expiration Date	12 months from date of receipt.



25 ug of the indicated Mouse tissue extr...



Lanes: Lane 1: 22 ug mouse liver cell ly...



WB Suggested Anti-PLD2 Antibody Titratio...