


Product Datasheet


DBI antibody (orb574529)

Description	Rabbit polyclonal antibody to DBI
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human DBI
Target	DBI
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MSQAEFEKAAEEVRHLKTKPSDEEMLFYIGHYKQATVGDINTERPGMLDF
Clonality	Polyclonal
MW	10kDa
Uniprot ID	P07108
NCBI	NP_001073331
Expiration Date	12 months from date of receipt.




Host: Rabbit, Target Name: DBI, Sample T...

DBI



Host: Rabbit, Target Name: DBI, Sample T...



Immunohistochemistry with Brain, cerebel...