

Product Datasheet

ZBTB4 antibody (orb574550)

Description	Rabbit polyclonal antibody to ZBTB4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human ZBTB4
Target	ZBTB4
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RPGGTPTPVIAYSKGSAGTRPGDVKEEAPQEMQVSSSSGEAGGGSTAAEE
Clonality	Polyclonal
MW	105kDa
Uniprot ID	Q9P1Z0
NCBI	NP_065950
Expiration Date	12 months from date of receipt.



WB Suggested
Anti-ZBTB4
Antibody
Titrati...