



Product Datasheet

Tead4 antibody (orb574665)

Description	Rabbit polyclonal antibody to Tead4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	Tead4
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RRKAREIQAKLKDQAAKNKALQSMAMSSAQIVSATAFHSKMALARGPGY
Clonality	Polyclonal
MW	48kDa
Uniprot ID	Q62296
NCBI	NP_035697
Expiration Date	12 months from date of receipt.



WB Suggested
Anti-Tead4
antibody
Titrati...



WB Suggested
Anti-Tead4
antibody
Titrati...



WB Suggested
Anti-Tead4
Antibody,
Titrati...