



Product Datasheet

PTHLH antibody (orb574866)


Description	Rabbit polyclonal antibody to PTHLH
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Goat, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human PTHLH
Target	PTHLH
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: YKEQLKTPGKKKKGKPGKRKEQEKKKRRTRSAWLDSGVTGSGLEGDHLS
Clonality	Polyclonal
MW	20kDa
Uniprot ID	P12272
NCBI	NP_002811
Expiration Date	12 months from date of receipt.



IHC
Information:Rabbit
Anti-PTHLH
Antibo...



Primary antibody
dilution: 1:1000,
Secon...



Rabbit Anti-PTHLH
Antibody, Paraffin
Emb...