


Product Datasheet


LIG4 antibody (orb574941)

Description	Rabbit polyclonal antibody to LIG4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human LIG4
Target	LIG4
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DGERMQMHKDGDVYKYFSRNGYNYTDQFGASPTGSLTPFIHNAFKADIQ
Clonality	Polyclonal
MW	104kDa
Uniprot ID	Q8IY66
NCBI	NP_002303
Expiration Date	12 months from date of receipt.

LIG4



Host: Rabbit,
Target Name:
LIG4, Sample ...



Human
Pancreas



WB Suggested
Anti-LIG4
Antibody
Titration...