


Product Datasheet

HOXD9 antibody (orb575086)

Description	Rabbit polyclonal antibody to HOXD9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human HOXD9
Target	HOXD9
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RGPSPGPGSPANGRHYGIKPETRAAPAPATAASTTSSSSTLSSSSKRTE
Clonality	Polyclonal
MW	35kDa
Uniprot ID	P28356
NCBI	NP_055028
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target: HOXD9,
Positive co...



WB Suggested
Anti-HOXD9
Antibody
Titration...