



Product Datasheet

YBX2 antibody (orb575091)

Description	Rabbit polyclonal antibody to YBX2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human YBX2
Target	YBX2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RKSRRFIPRPPSVAPPPMVAEIPSAGTGP GSKGERAEDSGQRPRRWCPPP
Clonality	Polyclonal
MW	39kDa
Uniprot ID	Q9Y2T7
NCBI	NP_057066
Expiration Date	12 months from date of receipt.



Human kidney



WB Suggested
Anti-YBX2
Antibody
Titratio...