



Product Datasheet

PSIP1 antibody (orb575182)


Description	Rabbit polyclonal antibody to PSIP1
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Guinea pig, Human, Mouse, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human PSIP1
Target	PSIP1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: AADRRKRKQEEQMETEQNKDEGKKPEVKKVEKKRETSMD SRLQRIHAEIK
Clonality	Polyclonal
MW	60kDa
Uniprot ID	O75475
NCBI	NP_150091
Expiration Date	12 months from date of receipt.



Human Jurkat



Rabbit Anti-PSIP1 antibody, Catalog Numb...



WB Suggested Anti-PSIP1 Antibody, Titrat...