


Product Datasheet

KCNA10 antibody (orb575397)

Description	Rabbit polyclonal antibody to KCNA10
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human KCNA10
Target	KCNA10
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DVCGWKEMEVALVNFDNSDEIQEEPGYATDFDSTSPKGRPGGSSFSNGKI
Clonality	Polyclonal
MW	58kDa
Uniprot ID	Q16322
NCBI	NP_005540
Expiration Date	12 months from date of receipt.



Anti-KCNA10 antibody IHC of human heart....



WB Suggested Anti-KCNA10 Antibody Titrat...