



Product Datasheet

TRPM4 antibody (orb575427)


Description	Rabbit polyclonal antibody to TRPM4
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human TRPM4
Target	TRPM4
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ELLTVYSSDGDSEEFETIVLKALVKACGSSEASAYLDELRLAVAWNRVDI
Clonality	Polyclonal
MW	134 kDa
Uniprot ID	Q8TD43
NCBI	NP_060106
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



Host: Rabbit,
Target Name:
TRPM4,
Sample...



Host: Rabbit,
Target Name:
TRPM4,
Sample...