


Product Datasheet


TPCN1 antibody (orb575428)

Description	Rabbit polyclonal antibody to TPCN1
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human TPCN1
Target	TPCN1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: YQEAAIYLQEGENNDKFFTHPKDAKALAAAYLFAHNHLFYLMELATALLLL
Clonality	Polyclonal
MW	94 kDa
Uniprot ID	Q9ULQ1
NCBI	NP_060371
Expiration Date	12 months from date of receipt.




25 ug of the indicated Human whole cell ...

TPCN1



Host: Rabbit,
Target Name:
TPCN1,
Sample...



Host: Rabbit,
Target Name:
TPCN1,
Sample...