


Product Datasheet

NMUR2 antibody (orb575438)

Description	Rabbit polyclonal antibody to NMUR2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human NMUR2
Target	NMUR2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MSGMEKLQNASWIYQKLEDPFQKHLNSTEEYLAFLCGPRRSHFFLPVSV
Clonality	Polyclonal
MW	46kDa
Uniprot ID	Q9GZQ4
NCBI	NP_064552
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target: NMUR2,
Positive co...



WB Suggested
Anti-NMUR2
Antibody,
Titrat...