



Product Datasheet

P2RX7 antibody (orb575528)


Description	Rabbit polyclonal antibody to P2RX7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human P2RX7
Target	P2RX7
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: IQSMNYGTIKWFFHVIIFSVCALVSDKLYQRKEPVISSVHTKVKGIAE
Clonality	Polyclonal
MW	69kDa
Uniprot ID	Q99572
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
P2RX7,
Sample...



Host: Rabbit,
Target: P2RX7,
Positive co...



Human kidney