



Product Datasheet

Atp6v0a1 antibody (orb576251)


Description	Rabbit polyclonal antibody to Atp6v0a1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide corresponding to a region of Mouse
Target	Atp6v0a1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TLNFGGIRVGNPTEEDAEIIQHDQLSTHSEDAEPEPTDEVDFDGMTMVH
Clonality	Polyclonal
MW	96 kDa
Uniprot ID	A2A5A0
NCBI	NP_058616
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



Application:
Western blotting,
Species+t...



Host: Rabbit,
Target Name:
Atp6v0a1,
Sam...