



Product Datasheet

HIF1A antibody (orb576507)

Description	Rabbit polyclonal antibody to HIF1A
Species/Host	Rabbit
Reactivity	Gallus, Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human HIF1A
Target	HIF1A
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VKGCKSSSEQNGMEQKTIILIPSDLACRLLGQSMDESGLPQLTSDYDCEVNA
Clonality	Polyclonal
MW	93kDa
Uniprot ID	Q16665
NCBI	NP_001521
Expiration Date	12 months from date of receipt.




25 ug of the indicated Mouse whole tissu...



HIF1A antibody - middle region (orb57650...)

Av-54 Av-54 BD
BSA Milk



HIF1A antibody - middle region (orb57650...)