


Product Datasheet

Egln2 antibody (orb577284)

Description	Rabbit polyclonal antibody to Egln2
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	Egln2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: FWSDRRNPEVKPAYATRYAITVWYFDAKERAARDKYQLASGQKGVQVP
Clonality	Polyclonal
MW	45kDa
Uniprot ID	Q6AYU4
NCBI	NP_001004083
Expiration Date	12 months from date of receipt.



Rabbit Anti-Egln2 antibody, Formalin Fix...



WB Suggested Anti-Egln2 Antibody, Titrat...