


Product Datasheet


RPL13 antibody (orb577492)

Description	Rabbit polyclonal antibody to RPL13
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human RPL13
Target	RPL13
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KKEKARVITEEEKNFKAFASLRMARANARLFGIRAKRAKEAAEQDVEKKK
Clonality	Polyclonal
MW	24kDa
Uniprot ID	P26373
NCBI	NP_000968
Expiration Date	12 months from date of receipt.

RPL13




Host: Rabbit,
Target Name:
RPL13,
Sample...



Host: Rabbit,
Target Name:
RPL13,
Sample...

RPL13



Host: Rabbit,
Target Name:
RPL13,
Sample...