



Product Datasheet

PUM3 antibody (orb577539)


Description	Rabbit polyclonal antibody to PUM3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human KIAA0020
Target	PUM3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GKKGVKQFKNKQGGDKSPKNKFQANKFNKKRKFQPDGRSDESAAKPKW
Clonality	Polyclonal
MW	71kDa
Uniprot ID	Q15397
NCBI	NP_001026861
Expiration Date	12 months from date of receipt.



Paraffin Embedded Tissue: Human Lung Cel...



Rabbit Anti-KIAA0020 Antibody, Paraffin ...



WB Suggested Anti-KIAA0020 Antibody Titr...