



## Product Datasheet

### HNRNPUL1 antibody (orb577738)


<b>Description</b>	Rabbit polyclonal antibody to HNRNPUL1
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, IP, WB
<b>Immunogen</b>	The immunogen is a synthetic peptide directed towards the C terminal region of human HNRNPUL1
<b>Target</b>	HNRNPUL1
<b>Preservatives</b>	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
<b>Form/Appearance</b>	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Synthetic peptide located within the following region: TYPQPSYNQYQQYAQQWNQYYQNQGQWPPYYGNIDYGSYSGNTQGGTSTQ
<b>Clonality</b>	Polyclonal
<b>MW</b>	83kDa
<b>Uniprot ID</b>	<a href="#">Q9BUJ2</a>
<b>NCBI</b>	<a href="#">NP_653333</a>
<b>Expiration Date</b>	12 months from date of receipt.



Amount and Sample Type:  
Lane 1: 5% Input...



Lane 1: 5 ug HeLa S3, Lane 2: 10 ug HeLa...



Rabbit Anti-HNRNPUL1 Antibody, Paraffin ...