


Product Datasheet

Wnt11 antibody (orb577976)

Description	Rabbit polyclonal antibody to Wnt11
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of Mouse Wnt11
Target	Wnt11
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PGCSCGPVPGPEPPGPNRWGGCADNLSYGLLMGAKFSDAPMKVKKTGSQA
Clonality	Polyclonal
MW	39kDa
Uniprot ID	P48615
NCBI	NP_033545
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
Wnt11,
Sample...



Rabbit Anti-
Wnt11 antibody,
Formalin Fix...