



Product Datasheet

WNT7B antibody (orb577994)


Description	Rabbit polyclonal antibody to WNT7B
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human WNT7B
Target	WNT7B
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: WTTLPKFREVGHLLKEKYNAAVQVEVVRASRLRQPTFLRIKQLRSYQKPM
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	P56706
NCBI	NP_478679
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



Anti-WNT7B antibody IHC of human brain, ...



Host: Mouse,
Target Name:
WNT7B, Sample
...