



Product Datasheet

LSS antibody (orb578239)

Description	Rabbit polyclonal antibody to LSS
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human LSS
Target	LSS
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TEGTCLRRRGGPYKTEPATDLGRWRLNCERGRQTWTYLQDERAGREQTGL
Clonality	Polyclonal
MW	83 kDa
Uniprot ID	P48449
NCBI	NP_002331
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



LSS antibody - N-terminal region (orb578...



WB Suggested Anti-LSS Antibody Titration...