



Product Datasheet

ALPP antibody (orb578290)

Description	Rabbit polyclonal antibody to ALPP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human ALPP
Target	ALPP
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TVLLYGNGPGYVLKDGARPDVTESESGSPEYRQQSAVPLDEETHAGEDVA
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	P05187
NCBI	NP_001623
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



Rabbit Anti-ALPP Antibody, Paraffin Embe...



WB Suggested Anti-ALPP Antibody Titratio...