


Product Datasheet

PLD3 antibody (orb579283)

Description	Rabbit polyclonal antibody to PLD3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human PLD3
Target	PLD3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: WVLLVLLILAVVGFALMTQLFLWEYGDHLHFGPNQRPA PCYDPCEAVLVE
Clonality	Polyclonal
MW	55 kDa
Uniprot ID	Q8IV08
NCBI	NP_001026866
Expiration Date	12 months from date of receipt.


PLD3



Host: Rabbit,
Target Name:
PLD3, Sample


...

PLD3



Host: Rabbit,
Target Name:
PLD3, Sample

...



Host: Rabbit,
Target: PLD3,
Positive con...