



Product Datasheet

ATP6V0C antibody (orb579396)


Description	Rabbit polyclonal antibody to ATP6V0C
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human ATP6V0C
Target	ATP6V0C
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VVAVLIANSLNDDISLYKSFLQLGAGLSVGLSGLAAGFAIGIVGDAGVRG
Clonality	Polyclonal
MW	16 kDa
Uniprot ID	P27449
NCBI	NP_001685
Expiration Date	12 months from date of receipt.



25 ug of the indicated whole cell extrac...



Host: Rabbit,
Target Name:
ATP6V0C,
Samp...



Host: Rabbit,
Target Name:
ATP6V0C,
Samp...