

Product Datasheet

Cdh7 antibody (orb579787)

Description	Rabbit polyclonal antibody to Cdh7
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Rat Cdh7
Target	Cdh7
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: WNQFFVLEEYMGSDPLYVGKLSHSDVDKGDGFIKYILSGEGASSIFIIDEN
Clonality	Polyclonal
MW	86kDa
Uniprot ID	Q5DWW2
NCBI	NP_001012755
Expiration Date	12 months from date of receipt.



WB Suggested
Anti-Cdh7
Antibody,
Titrati...