



Product Datasheet

COX15 antibody (orb579790)


Description	Rabbit polyclonal antibody to COX15
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human COX15
Target	COX15
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DWHLIKEMKPPTSQEEWEAEFQRYQQFPEFKILNHDMTLTFKFIWYMEY
Clonality	Polyclonal
MW	46 kDa
Uniprot ID	Q7KZN9
NCBI	NP_004367
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



Lanes: Lane 1: 20 ug human fibroblast mi...



Rabbit Anti-COX15 Antibody, Paraffin Emb...