



Product Datasheet

Tmed1 antibody (orb579880)

Description	Rabbit polyclonal antibody to Tmed1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide corresponding to a region of Rat
Target	Tmed1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EAGPPPIQDGEFTFLLPAGRKQCFYQSAPANASLETEYQVIGGAGLDVDF
Clonality	Polyclonal
MW	25kDa
Uniprot ID	Q5BK85
NCBI	NP_001013450
Expiration Date	12 months from date of receipt.



Lanes: Lane 1 and 2: 30 ug HEK-293 cell ...



Rabbit Anti-Tmed1 Antibody, Catalog Numb...



WB Suggested Anti-Tmed1 antibody Titrati...