



Product Datasheet

EEF1A1 antibody (orb580310)

Description	Rabbit polyclonal antibody to EEF1A1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human EEF1A1
Target	EEF1A1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: IVDMVPGKPMCVESFSDYPLGRFAVRDMRQTVAVGVKAVDKKAAGAGK
Clonality	Polyclonal
MW	50kDa
Uniprot ID	P68104
NCBI	NP_001393
Expiration Date	12 months from date of receipt.



Host: Mouse,
Target Name:
EEF1A1,
Sample...



WB Suggested
Anti-EEF1A1
antibody
Titrat...



WB Suggested
Anti-EEF1A1
antibody
Titrat...