



Product Datasheet

PFN1 antibody (orb580352)


Description	Rabbit polyclonal antibody to PFN1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human PFN1
Target	PFN1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: AGWNAYIDNLMADGTCQDAIVGYKDSPSVWAAVPGKTFVNITPAEVGVL
Clonality	Polyclonal
MW	15 kDa
Uniprot ID	P07737
NCBI	NP_005013
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



Lanes: 1. 20 ug human pregnant uterine m...



Lanes: 1. 20 ug murine non-pregnant uter...