



Product Datasheet

PNPLA3 antibody (orb580585)

Description	Rabbit polyclonal antibody to PNPLA3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human PNPLA3
Target	PNPLA3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: CSPKGCPAETKAEATPRSLRSSLNFFLGKVPAGAEGLSTFPSFSLEKS
Clonality	Polyclonal
MW	53kDa
Uniprot ID	Q9NST1
NCBI	NP_079501
Expiration Date	12 months from date of receipt.




Host: Rabbit,
Target Name:
PNPLA3,
Sampl...



Host: Rabbit,
Target Name:
PNPLA3,
Sampl...

PNPLA3



Host: Rabbit,
Target Name:
PNPLA3,
Sampl...