



Product Datasheet

UBA3 antibody (orb581104)


Description	Rabbit polyclonal antibody to UBA3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human UBA3
Target	UBA3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EGRWNHVKKFLERSGPFTHPDFEPSTESLQFLDTCCKVLVIGAGGLGCEL
Clonality	Polyclonal
MW	52kDa
Uniprot ID	Q8TBC4
NCBI	NP_003959
Expiration Date	12 months from date of receipt.



Sample Type: 1.
Human NT-2
cells (60 ug)...



Sample Type:
Human brain
stem cells, Pri...



WB Suggested
Anti-UBA3
Antibody
Titratio...