



Product Datasheet

Ift88 antibody (orb582172)

Description	Rabbit polyclonal antibody to Ift88
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Target	Ift88
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MMENVHLAPETDEDDLYSGFNDYNPAYDTELENDTGFQQAVRTSHGRRP
Clonality	Polyclonal
MW	89kDa
Uniprot ID	D4ACI9
NCBI	157821925 , B5WUVB49013 , NP_001100736
Expiration Date	12 months from date of receipt.



Host: Rabbit, Target Name: IFT88, Sample...



Immunohistochemistry with Testis tissue.



WB Suggested Anti-Ift88 Antibody Titration...