



Product Datasheet

GNAI2 antibody (orb582416)


Description	Rabbit polyclonal antibody to GNAI2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human GNAI2
Target	GNAI2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EYTGANKYDEAASYIQSKFEDLNKRKDTKEIYTHFTCATDTKNVQFVFDA
Clonality	Polyclonal
MW	40kDa
Uniprot ID	P04899
NCBI	NP_002061
Expiration Date	12 months from date of receipt.



Sample Type:
Nthy-ori cell
lysate (50 ug...



Host: Mouse,
Target Name:
GNAI2, Sample
...



WB Suggested
Anti-GNAI2
Antibody
Titration...