



Product Datasheet

UBQLN1 antibody (orb582568)


Description	Rabbit polyclonal antibody to UBQLN1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human UBQLN1
Target	UBQLN1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: QFGGNPFASLVSNNTSSGEGSQPSRTENRDPLPNPWAPQTSQSSSASSGTA
Clonality	Polyclonal
MW	62kDa
Uniprot ID	Q5T6J9
NCBI	NP_038466
Expiration Date	12 months from date of receipt.



Anti-UBQLN1 antibody IHC of human brain,...



Anti-UBQLN1 antibody IHC of human prosta...



Immunohistochemistry with Human Lung, re...