

Product Datasheet

Sec14I3 antibody (orb582859)

Description	Rabbit polyclonal antibody to Sec14I3
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide corresponding to a region of Rat
Target	Sec14I3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RDQVKTYEHSVQISRGSSHQVEYEILFPGCVLRWQFSSDGADIGFGVFL
Clonality	Polyclonal
MW	44kDa
Uniprot ID	Q9Z1J8
NCBI	NP_072130
Expiration Date	12 months from date of receipt.



WB Suggested
Anti-Sec14I3
Antibody, Titr...