



Product Datasheet

UCHL5 antibody (orb583318)


Description	Rabbit polyclonal antibody to UCHL5
Species/Host	Rabbit
Reactivity	Human, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human UCHL5
Target	UCHL5
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DGLREGPIDLGACNQDDWISAVRPVIEKRIQKYSEGEIRFNLMAIVSDRK
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	Q9Y5K5
NCBI	NP_057068
Expiration Date	12 months from date of receipt.



25 ug of the indicated Human whole cell ...



WB Suggested Anti-UCHL5 Antibody, Positi...



WB Suggested Anti-UCHL5 Antibody Titration...