


Product Datasheet

DEFB129 antibody (orb583999)


Description	Rabbit polyclonal antibody to DEFB129
Species/Host	Rabbit
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human DEFB129
Target	DEFB129
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TSFFANTNFVIIPNATPMNSATISTMTPGQITYTATSTKSNTKESRDSAT
Clonality	Polyclonal
MW	18kDa
Uniprot ID	Q9H1M3
NCBI	NP_543021
Expiration Date	12 months from date of receipt.



Host: Rabbit,
Target Name:
DEFB129,
Samp...



Rabbit Anti-
DEFB129
antibody,
Formalin F...



Sample Type:
Pig bronchus
smooth muscle

...