

Product Datasheet

LRRFIP2 antibody (orb584199)

Description	Rabbit polyclonal antibody to LRRFIP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	LRRFIP2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RKLKLQLEEERQKCSRNDGTVGDLAGLQNGSDLQFIEMQRDANRQISEYK
Clonality	Polyclonal
MW	82kDa
Uniprot ID	Q9Y608
NCBI	NP_006300
Expiration Date	12 months from date of receipt.



WB Suggested
Anti-LRRFIP2
Antibody, Titr...