



Product Datasheet

USP2 antibody (orb584295)


Description	Rabbit polyclonal antibody to USP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human USP2
Target	USP2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: NEVNRVTLRPKSNPENLDHLPDDEKGRQMWRKYLEREDSRIGDLFVGQLK
Clonality	Polyclonal
MW	68kDa
Uniprot ID	O75604
NCBI	NP_004196
Expiration Date	12 months from date of receipt.



Rabbit Anti-USP2 antibody, Formalin Fixe...



Rabbit Anti-USP2 antibody, Formalin Fixe...



Rabbit Anti-USP2 Antibody, Catalog Numbe...