



Product Datasheet

TAAR5 antibody (orb584397)


Description	Rabbit polyclonal antibody to TAAR5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human TAAR5
Target	TAAR5
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TTLSKSLAGAAKHERKAAKTLGIAVGIYLLCWLPFTIDTMVDSLLHFITP
Clonality	Polyclonal
MW	38kDa
Uniprot ID	O14804
NCBI	NP_003958
Expiration Date	12 months from date of receipt.



Immunofluorescence
- dilution: 1.3
ug/ml...



WB Suggested Anti-
TAAR5 Antibody
Titration...



WB Suggested Anti-
TAAR5 antibody
Titration...