

Product Datasheet

PBK/TOPK (phospho-Thr9) antibody (orb5844)

Description	PBK/TOPK (phospho-Thr9) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PBK/TOPK around the phosphorylation site of Thr9 (FK(p-T)P)
Target	Phospho-PBK (Thr9)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	Q96KB5
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.